



EON Inspiring for Splendid Healthcare®
MEDITECH
PVT. LTD.

CATALOGUE



Otorhinolaryngology Devices Par Excellence

CDSCO
Central Drugs Standard
Control Organization

ISO
9001:2015

ISO
13485:2016

CE

SFDA

ANVISA

TGA
CERTIFIED

UK MHRA
Regulating Medicines & Medical

US FDA
U.S. Food and Drug Administration



Tita - Prostheses type Partial Vario Campana

(Replacement of Malleus & Incus)

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
E1101	0.30	3.00	5.00-1.75	5.00	1 Unit



Tita - Prostheses type Total Vario Cask

(Replacement of Complete Ossicular chain)

Material : CP Titanium Medical Grade ASTM F 67


REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
E1102	0.30	3.00	7.00-3.00	7.00	1 Unit



Titanium Universal Middle Ear Prosthesis

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic

REF Number	Shaft Diameter	Disk Diameter	Connector Diameter	Total Length	Units/Box Gamma Sterile
E1138	0.30	3.00	1.00	7.40	1 Unit



TitaHap prosthesis type PORP & TORP

Material : CP Titanium Medical Grade ASTM F 67 & HAp

REF Number	Description	Units/Box Gamma Sterile
E1207	TitaHap prosthesis type Vario Partial (PORP)	1 Unit
E1208	TitaHap prosthesis type Vario Total (TORP)	1 Unit



Clip Stapes Piston

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
E1121	7.00 to 3.50	0.40	1 Unit



Tita - prosthesis Piston type Loop

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
E1103	7.00 to 3.50	0.40	1 Unit



Tita - PTFE Vario Piston type Loop

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic(PTFE)

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
E4119	7.00 to 4.00	0.40	1 Unit

All dimensions are in mm



Super Titanium Vario Piston type Loop

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
E9131	7.00 to 3.50	0.40	1 Unit



Super Tita - PTFE Vario Piston type Loop

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic (PTFE)

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
E9132	7.00 to 4.00	0.40	1 Unit



Super Elastic NiTiNOL Vario Piston type Loop

Material : NiTiNOL & Fluoroplastic

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
E5120	7.00 to 3.50	0.40	1 Unit



Super Elastic NiTiNOL - PTFE Vario Piston type Loop

Material : NiTiNOL & Fluoroplastic (PTFE)

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
E8130	7.00 to 4.00	0.40	1 Unit



Titanium Stapes Piston with Integrated PTFE Shoe

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic (PTFE)

REF Number	Functional Length	Total Length	Units/Box Gamma Sterile
E1122	7.00 to 3.50	7.00	1 Unit



Fluoroplastic Partial Prosthesis (PORP)

Material : Fluoroplastic

REF Number	Shaft Diameter OD/ID	Disk Diameter	Functional Length	Total Length	Units/Box ETO Sterile
E2133	1.8/1.20	3.00	5.00-2.00	5.00	1 Unit



Fluoroplastic Total Prosthesis (TORP)

Material : Fluoroplastic

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box ETO Sterile
E2134	0.8	3.00	7.00-3.00	7.00	1 Unit



Fluoroplastic Offset Partial Prosthesis (PORP)

Material : Fluoroplastic

REF Number	Shaft Diameter OD/ID	Functional Length	Total Length	Units/Box ETO Sterile
E2141	1.8/1.20	5.00-2.00	5.00	1 Unit



Fluoroplastic Offset Total Prosthesis (TORP)

Material : Fluoroplastic

REF Number	Shaft Diameter	Functional Length	Total Length	Units/Box ETO Sterile
E2142	0.80	7.00-3.00	7.00	1 Unit



Fluoroplastic Universal Middle Ear Prosthesis Offset

Material : Fluoroplastic

In case of TORP-Shaft OD 0.8mm | In case of PORP-Shaft ID 1.20mm

REF Number	Functional Length	Units/Box ETO Sterile
E2139	7.00 to 2.00	1 Unit



Fluoroplastic Universal Middle Ear Prosthesis Centered

Material : Fluoroplastic

In case of TORP-Shaft OD 0.8mm | In case of PORP-Shaft ID 1.20mm

REF Number	Functional Length	Units/Box ETO Sterile
E2140	7.00 to 2.00	1 Unit



Fluoroplastic Stapes Piston

Material : Fluoroplastic

REF Number	Shaft Diameter	Loop Inner Diameter	Functional Length	Total Length	Units/Box ETO Sterile
E2127	0.5	0.60	6.00-3.50	6.00	1 Unit
E2128	0.4	0.60	6.00-3.50	6.00	1 Unit
E2129	0.6	0.60	6.00-3.50	6.00	1 Unit



Middle Ear Implant Cutter

Material : Stainless Steel

REF Number	Description	Units/Box
EI01	Middle Ear Implant Cutter	1 Unit, Non-Sterile



Middle Ear Implant Sizer

Material : Stainless Steel

REF Number	Description	Units/Box Non-Sterile
EI02	Middle Ear Implant Sizer	1 Unit

All dimensions are in mm



Tita - prosthesis vent tube type collar button

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
E1104	1.40	2.55	2.55	1.20	10 Units
E1105	1.30	2.55	2.55	1.20	10 Units
E1106	1.50	2.75	2.75	1.20	10 Units

Shepard Titanium

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
E1151	1.14	2.40	2.40	2.20	10 Units

Armstrong | Armstrong Blue | Armstrong with Tail

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2109	1.14	3.54	2.40	3.85	10 Units
E2109B	1.14	3.54	2.40	3.85	10 Units
E2109T	1.14	3.54	2.40	3.85	10 Units

Modified Armstrong | Modified Armstrong Blue Modified Armstrong with Tail

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2110	1.14	3.46	2.50	3.04	10 Units
E2110B	1.14	3.46	2.50	3.04	10 Units
E2110T	1.14	3.46	2.50	3.04	10 Units

Donaldson

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2124	1.14	2.3	2.3	2.2	10 Units
E2124B	1.14	2.3	2.3	2.2	10 Units

Reuter Bobbin

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2137	1.14	2.3	2.3	1.6	10 Units
E2137B	1.14	2.3	2.3	1.6	10 Units

All dimensions are in mm



Shah and Mini Shah

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2112	1.14	3.4	2.3	2.4	10 Units
E2112B	1.14	3.4	2.3	2.4	10 Units
E2112T	1.14	3.4	2.3	2.4	10 Units
E2112W	1.14	3.4	2.3	2.4	10 Units
E2112BW	1.14	3.4	2.3	2.4	10 Units
E2112B	1.14	3.4	2.3	2.4	10 Units
E2113	0.76	2.2	1.5	1.5	10 Units
E2113B	0.76	2.2	1.5	1.5	10 Units
E2113T	0.76	2.2	1.5	1.5	10 Units
E2113W	0.76	2.2	1.5	1.5	10 Units

Shepard & Mini Shepard

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2114	1.14	2.4	2.4	2.2	10 Units
E2114B	1.14	2.4	2.4	2.2	10 Units
E2114T	1.14	2.4	2.4	2.2	10 Units
E2114W	1.14	2.4	2.4	2.2	10 Units
E2114BW	1.14	2.4	2.4	2.2	10 Units
E2114BT	1.14	2.4	2.4	2.2	10 Units
E2138	1.00	2.4	2.4	2.2	10 Units
E2138W	1.00	2.4	2.4	2.2	10 Units
E2138T	1.00	2.4	2.4	2.2	10 Units
E2138B	1.00	2.4	2.4	2.2	10 Units
E2138BW	1.00	2.4	2.4	2.2	10 Units
E2138BT	1.00	2.4	2.4	2.2	10 Units

Collar Button (ID 1.27 mm) & (ID 1.14 mm)

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2123	1.27	3.0	3.0	2.15	10 Units
E2123B	1.27	3.0	3.0	2.15	10 Units
E2123T	1.27	3.0	3.0	2.15	10 Units
E2125	1.14	2.50	2.50	2.0	10 Units
E2125B	1.14	2.50	2.50	2.0	10 Units
E2125T	1.14	2.50	2.50	2.0	10 Units

All dimensions are in mm



Baxter Beveled Bobbin

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2136	1.27	2.8	2.8	2.2	10 Units
E2136B	1.27	2.8	2.8	2.2	10 Units



Bevel Bobbin

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2117	1.14	2.8	2.8	1.7	10 Units
E2117B	1.14	2.8	2.8	1.7	10 Units



Straight Tube

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2116	1.14	2.7	1.54	7.25	10 Units



Donaldson Silicone

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E3111	1.10	2.3	2.3	2.2	10 Units



T - Tube

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Total Length	Units/Box ETO Sterile
E3115S	1.14	8.0	6.0	10 Units
E3115M	1.14	8.0	9.0	10 Units
E3115L	1.14	8.0	12.0	10 Units



Paparella

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E3126	1.10	3.5	2.5	1.65	10 Units



Gold vent Tube Type Collar Button With Wire

Material : Gold Plated Silver


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
E1170	1.25	2.55	2.55	1.20	10 Units

All dimensions are in mm



Nasal Septal Button

Material : Medical grade Silicone



REF Number	Description	Size(cm)	Units/Box ETO Sterile
E3145S	Nasal Septal Button - Standard	3.0	1 Unit
E3145L	Nasal Septal Button - Large	5.0	1 Unit
E3145E	Nasal Septal Button - Eliptical	3.0	1 Unit
E3145EL	Nasal Septal Button - Eliptical Large	5.0	1 Unit

Internal Nasal Airway Splint


Material : Medical grade Silicone



REF Number	Descriptions	No. of ETO Sterile Pairs/Packages
E3147S	Internal Nasal Airway Splint - Standard	10 Pairs
E3147L	Internal Nasal Airway Splint - Large	10 Pairs

Internal Nasal Pre-Cut Splint

Material : Medical grade Silicone



REF Number	Descriptions	No. of ETO Sterile Pairs/Packages
E3147P	Internal Nasal Pre-cut Splint	10 Pairs

Internal Bi-Valve Splints

Material : Medical grade Silicone & Fluoroplastic



REF Number	Descriptions	Thickness (mm)	Units/Box ETO Sterile
E3148SS	Silicone-standard	0.50	10 Pairs
E3148ST	Silicone-Thick	1.00	10 Pairs
E3148FS	Fluoroplastic-Standard	0.50	10 Pairs
E3148FT	Fluoroplastic-Thick	0.25	10 Pairs
E3148FSP	Pediatric-Standard	0.50	10 Pairs
E3148FTP	Pediatric-Thin	0.25	10 Pairs

External Thermoplastic Splints



REF Number	Descriptions	Units/Box ETO Sterile/Non Sterile
E5150S	External Thermoplastic Splint Small	10 Units
E5150M	External Thermoplastic Splint Medium	10 Units
E5150L	External Thermoplastic Splint Large	10 Units

External Aluminium Splints



REF Number	Descriptions	Units/Box ETO Sterile/Non Sterile
E6149S	External Aluminium Splints Small	10 Units
E6149M	External Aluminium Splints Medium	10 Units
E6149L	External Aluminium Splints Large	10 Units



Silicone Sheets

Material : Medical grade Silicone



REF Number	Size (mm)	Thickness (mm)	Units/Box ETO Sterile
E3146S	55x75	3.00	10 Units
E3146T	55x75	5.00	10 Units
E3146NS	55x75	1.00	10 Units
E3146MS	55x75	0.50	10 Units
E3146US	55x75	0.25	10 Units
E3146UUS	38x50	0.13	10 Units

Silicone Blocks

Material: Medical Grade Silicone



REF Number	Size (mm)	Thickness (mm)	Units/Box ETO Sterile
E3146BS	50x30x10	10.00	1 Unit
E3146BM	60x40x15	15.00	1 Unit
E3146BL	70x50x20	20.00	1 Unit

Epistaxis Balloon Catheter

Material: Medical Grade Silicone



REF Number	Product Name	Specification	Units/Box ETO Sterile
E31030	Intra Pack	30 CC / ml	10 Units
E31010	Post pack	10 CC / ml	10 Units
E31031	Double Balloon	30 & 10 CC / ml	10 Units

External Splint Kits

Aluminium Splint kit

E6149K



Aluminium Splint (Small) OR	1 Unit
Aluminium Splint (Medium) OR	1 Unit
Aluminium Splint (Large)	1 Unit
Alcohol Swab	1 Unit
Skin Preparation Swab	1 Unit
Skin Closure Tape	1 Unit
Dorsum Pad	1 Unit

Thermoplastic Splint kit

E5150K



Thermoplastic Splint (Small) OR	1 Unit
Thermoplastic Splint (Medium) OR	1 Unit
Thermoplastic Splint (Large)	1 Unit
Alcohol Swab	1 Unit
Skin Preparation Swab	1 Unit
Skin Closure Tape	1 Unit
Dorsum Pad	1 Unit



Ear Wick

REF Number	Product Name	Specification	Product Description	Units/Box Gamma Sterile
E7161	Ear Wick	Dia.9.0mmX15.0mm	No String	20 Units
E7162	Ear Wick	Dia.9.0mmX24.0mm	No String	20 Units
E7163	Ear Pack	Dia.12.0mmX15.0mm	No String	20 Units
E7164	Ear Pack	Dia.12.0mmX24.0mm	No String	20 Units
E7169	Ear Wick	Dia.9.0mmX15.0mm	Fenestrated	20 Units
E7170	Ear Wick	Dia.9.0mmX24.0mm	Fenestrated	20 Units

Ear Packs

REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7165	Dia.9.0mmX15.0mm	With String	20 Units
E7166	Dia.9.0mmX24.0mm	With String	20 Units
E7167	Dia.12.0mmX15.0mm	With String	20 Units
E7168	Dia.12.0mmX24.0mm	With String	20 Units

Anatomical Nasal Dressings

REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7188	8.0cmX1.5cmX3.0cm	With String	10 Units
E7189	8.0cmX1.5cmX3.0cm	With String & Tube	10 Units
E7190	8.0cmX1.0cmX2.0cm	Thin	10 Units

Sinus packs

REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7191	3.5cmX1.2cmX0.6cm	With String	20 Units
E7192	3.5cmX1.2cmX0.9cm	With String	20 Units


Standard Nasal Dressings

REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7176	4.5cmX1.5cmx2.0cm	No String	10 Units
E7177	4.5cmX1.5cmx2.0cm	With String	10 Units
E7178	4.5cmX1.5cmx2.0cm	With String & Tube	10 Units
E7179	4.5cmX1.0cmx2.0cm	Thin	10 Units
E7180	8.0cmX1.5cmx2.0cm	No String	20 Units
E7181	8.0cmX1.5cmx2.0cm	With String	20 Units
E7182	8.0cmX1.5cmx2.0cm	With String & Tube	10 Units
E7183	8.0cmX1.0cmx2.0cm	Thin	10 Units



Epistaxis Nasal Dressings

REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7184	10.0cmX1.5cmX2.5cm	No String	10 Units
E7185	10.0cmX1.5cmX2.5cm	With String	10 Units
E7186	10.0cmX1.5cmX2.5cm	With String & Tube	10 Units
E7187	10.0cmX1.0cmX2.5cm	Thin	10 Units
E7193	5.5cmX1.0cmX2.5cm	Pediatric	10 Units



CoatAid With Proprietary Surface Treatment Technology


Advance Nasal Dressings

REF. Number	Specification	Product Description	Units/Box Gamma Sterile
C7181	8.0cmX1.5cmx2.0cm	With String	20 Units
C7185	10.0cmX1.5cmx2.5cm	With String	20 Units



BACILK Bacteriostatic Treatment Technology

REF. Number	Specification	Product Description	Units/Box Gamma Sterile
EB7181	8.0cmX1.5cmx2.0cm	With String	20 Units
EB7185	10.0cmX1.5cmx2.5cm	With String	20 Units



Lamina Laminated Structured Surface Technology

REF. Number	Specification	Product Description	Units/Box Gamma Sterile
L7181	8.0cmX1.5cmx2.0cm	With String	20 Units
L7185	10.0cmX1.5cmx2.5cm	With String	20 Units



MICROSPON Antimicrobial Treatment Technology

REF. Number	Specification	Product Description	Units/Box Gamma Sterile
EA7181	8.0cmX1.5cmx2.0cm	With String	20 Units
EA7185	10.0cmX1.5cmx2.5cm	With String	20 Units



Stapes Prosthesis Handling Forcep

Material : Stainless Steel

REF Number	Description	Units/Box Non-Sterile
EI03	Titanium-Stapes Prosthesis Handling Forceps	1 Unit
EI05	PTFE-Stapes prosthesis Handling Forceps	1 Unit



Stapes Prosthesis Crimping Forcep

Material : Stainless Steel

REF Number	Description	Units/Box Non-Sterile
EI04	Stapes Prosthesis Crimping Forceps	1 Unit



Tilley's Nasal Dressings Forcep

Material : Stainless Steel

REF Number	Description	Units/Box Non-Sterile
EI06	Tilley's Nasal Dressings Forceps	1 Unit



Micro Ear Alligator Forcep

Material : Stainless Steel

REF Number	Jaw Length	Length	Units/Box, Non-Sterile
EI0969	3.5mm	75mm	1 Unit



Ear Curettes

Material : Polyethylene


REF Number	Product Name	Size	Units/Box Non-Sterile
EI0962	Jobson Horne Probe-Round	3 mm	100 Units
EI0963	Jobson Horne Probe-Oval	2.5 mm	100 Units
EI0964	Cerumen Hook	Standard	100 Units



Suction Controller

Material : Polyethylene

REF Number	Length	Units/Box Non-Sterile
EI0968	39mm	10 Units




Disposable Suction Tube for Ear

Material : Stainless Steel & PE

10 Units/Box ETO Sterile

REF Number	Specification	REF Number	Specification
EI152580	Dia.2.5mmX/80.0mm	EI151080	Dia.1.0mmX80.0mm
EI153080	Dia.3.0mmX80.0mm	EI1510100	Dia.1.0mmX100.0mm
EI152060	Dia.2.0mmX60.0mm	EI150870	Dia.0.8mmX70.0mm
EI152070	Dia.2.0mmX70.0mm	EI1508100	Dia.0.8mmX100.0mm
EI152080	Dia.2.0mmX80.0mm	EI150760	Dia.0.7mmX60.0mm
EI1520100	Dia.2.0mmX100.0mm	EI150770	Dia.0.7mmX70.0mm
EI1516100	Dia.1.6mmX100.0mm	EI150780	Dia.0.7mmX80.0mm
EI151460	Dia.1.4mmX60.0mm	EI1507100	Dia.0.7mmX100.0mm
EI151470	Dia.1.4mmX70.0mm	EI150670	Dia.0.6mmX70.0mm
EI151480	Dia.1.4mmX80.0mm	EI150460	Dia.0.4mmX60.0mm
EI1514100	Dia.1.4mmX100.0mm	EI150470	Dia.0.4mmX70.0mm
EI151270	Dia.1.2mmX70.0mm	EI150480	Dia.0.4mmX80.0mm
EI1512100	Dia.1.2mmX100.0mm	EI1504100	Dia.0.4mmX100.0mm
EI151060	Dia.1.0mmX60.0mm	E10120	Ear Disposable Set
EI151070	Dia.1.0mmX70.0mm		



Reusable Suction Tube for Ear

Material : Stainless Steel

10 Units/Box ETO Sterile

REF Number	Specification	REF Number	Specification
EI140460	Dia.0.4mmX60.0mm	EI141470	Dia.1.4mmX70.0mm
EI140470	Dia.0.4mmX70.0mm	EI141480	Dia.1.4mmX80.0mm
EI140480	Dia.0.4mmX80.0mm	EI141670	Dia.1.6mmX70.0mm
EI140670	Dia.0.6mmX70.0mm	EI141770	Dia.1.7mmX70.0mm
EI140770	Dia.0.7mmX70.0mm	EI142070	Dia.2.0mmX70.0mm
EI140780	Dia.0.7mmX80.0mm	EI1420100	Dia.2.0mmX100.0mm
EI140870	Dia.0.8mmX70.0mm	EI142580	Dia.2.5mmX80.0mm
EI141070	Dia.1.0mmX70.0mm	EI143080	Dia.3.0mmX80.0mm
EI141270	Dia.1.2mmX70.0mm	EI141RS	Ear Reusable Set
EI1412100	Dia.1.2mmX100.0mm		



Disposable Suction Tube for Nose

Material : Stainless Steel & PE

10 Units/Box ETO Sterile


REF Number	Specification	REF Number	Specification
EI1220100	Dia.2.0mmX100.0mm	EI1235140	Dia.3.5mmX140.0mm
EI1220120	Dia.2.0mmX120.0mm	EI1235160	Dia.3.5mmX160.0mm
EI1220140	Dia.2.0mmX140.0mm	EI1240100	Dia.4.0mmX100.0mm
EI1220150	Dia.2.0mmX150.0mm	EI1240120	Dia.4.0mmX120.0mm
EI1220160	Dia.2.0mmX160.0mm	EI1240140	Dia.4.0mmX140.0mm
EI1230100	Dia.3.0mmX100.0mm	EI1240150	Dia.4.0mmX150.0mm
EI1230120	Dia.3.0mmX120.0mm	EI1240160	Dia.4.0mmX160.0mm
EI1230140	Dia.3.0mmX140.0mm	EI1250120	Dia.5.0mmX120.0mm
EI1230150	Dia.3.0mmX150.0mm	EI1250140	Dia.5.0mmX140.0mm
EI1230160	Dia.3.0mmX160.0mm	EI1250160	Dia.5.0mmX160.0mm
EI1235120	Dia.3.5mmX120.0mm	EI121DS	Nose Disposable Set



Reusable Suction Tube for Nose

Material : Stainless Steel


10 Units/Box ETO Sterile



REF Number	Specification	REF Number	Specification
EI1320120R	Dia.2.0mmX120.0mm	EI1335120L	Dia.3.5mmX120.0mm
EI1320120L	Dia.2.0mmX120.0mm	EI1335140R	Dia.3.5mmX140.0mm
EI1320140R	Dia.2.0mmX140.0mm	EI1335140L	Dia.3.5mmX140.0mm
EI1320140L	Dia.2.0mmX140.0mm	EI1335160R	Dia.3.5mmX160.0mm
EI1320160R	Dia.2.0mmX160.0mm	EI1335160L	Dia.3.5mmX160.0mm
EI1320160L	Dia.2.0mmX160.0mm	EI1340120R	Dia.4.0mmX120.0mm
EI1325120R	Dia.2.5mmX120.0mm	EI1340120L	Dia.4.0mmX120.0mm
EI1325120L	Dia.2.5mmX120.0mm	EI1340140R	Dia.4.0mmX140.0mm
EI1325140R	Dia.2.5mmX140.0mm	EI1340140L	Dia.4.0mmX140.0mm
EI1325140L	Dia.2.5mmX140.0mm	EI1340160R	Dia.4.0mmX160.0mm
EI1325160R	Dia.2.5mmX160.0mm	EI1340160L	Dia.4.0mmX160.0mm
EI1325160L	Dia.2.5mmX160.0mm	EI1350120R	Dia.5.0mmX120.0mm
EI1330120R	Dia.3.0mmX120.0mm	EI1350120L	Dia.5.0mmX120.0mm
EI1330120L	Dia.3.0mmX120.0mm	EI1350140R	Dia.5.0mmX140.0mm
EI1330140R	Dia.3.0mmX140.0mm	EI1350140L	Dia.5.0mmX140.0mm
EI1330140L	Dia.3.0mmX140.0mm	EI1350160R	Dia.5.0mmX160.0mm
EI1330160R	Dia.3.0mmX160.0mm	EI1350160L	Dia.5.0mmX160.0mm
EI1330160L	Dia.3.0mmX160.0mm	EI131RS	Nose Reusable Set
EI1335120R	Dia.3.5mmX120.0mm		

Disposable Throat Suction Tube


Material : Stainless Steel & PE



REF Number	Specification	Units/Box ETO Sterile
EI1620250	Dia.2.0mmX250.0mm	10 Units
EI1630250	Dia.3.0mmX250.0mm	10 Units

Disposable Sinus Suction Tube

Material : Stainless Steel & PE




REF Number	Specification	Units/Box ETO Sterile
EI1825150	Dia.2.5mmX150.0mm	10 Units
EI1830150	Dia.3.0mmX150.0mm	10 Units
EI1840150	Dia.4.0mmX150.0mm	10 Units

Ear Specula

Material : Polyethylene



REF Number	Tip Diameter	Units/Box Non-Sterile
EI0951	3mm	100 Units
EI0952	4mm	100 Units
EI0953	5mm	100 Units
EI0954	6mm	100 Units



Myringotomy Knife

Material : SS & Polyethylene

REF Number	Length	Units/Box, Non-Sterile
EI0966	157mm	1 Unit



Rosen Inserter

Material : Polyethylene


REF Number	Units/Box, Non-Sterile
EI0959	100 Units



ENT Burs

10 Units/Box, Non-Sterile

Stainless Steel	Tungsten Carbide	Fine Diamond	Course Diamond	Length (mm)	Head Diameter (mm)
EI2105SS	EI2105TC	EI2105FD	EI2105CD	70.0	0.5
EI2110SS	EI2110TC	EI2110FD	EI2110CD	70.0	1.0
EI2115SS	EI2115TC	EI2115FD	EI2115CD	70.0	1.5
EI2120SS	EI2120TC	EI2120FD	EI2120CD	70.0	2.0
EI2125SS	EI2125TC	EI2125FD	EI2125CD	70.0	2.5
EI2130SS	EI2130TC	EI2130FD	EI2130CD	70.0	3.0
EI2135SS	EI2135TC	EI2135FD	EI2135CD	70.0	3.5
EI2140SS	EI2140TC	EI2140FD	EI2140CD	70.0	4.0
EI2145SS	EI2145TC	EI2145FD	EI2145CD	70.0	4.5
EI2150SS	EI2150TC	EI2150FD	EI2150CD	70.0	5.0
EI2155SS	EI2155TC			70.0	5.5
EI2160SS	EI2160TC	EI2160FD	EI2160CD	70.0	6.0
EI2165SS	EI2165TC			70.0	6.5
EI2170SS	EI2170TC		EI2170CD	70.0	7.0



Straight Sinus Blades

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box, ETO Sterile/Non-Sterile
EI2201	4mm	0°	5000RPM	1 Unit



Curved Sinus Blades

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box, ETO Sterile/Non-Sterile
EI2202	4mm	40°	5000RPM	1 Unit
EI2203	4mm	60°	5000RPM	1 Unit



Tonsillectomy and Adenoidectomy Blades

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box, ETO Sterile/Non-Sterile
EI2204	4mm	40°	5000RPM	1 Unit



Myringotomy Kit



EI0972 Pediatric

Jobson Horne Probe
Ear Speculas- 3,4,5 MM
Myringotomy Knife
Suction Tube - 16,18 Gauge
Micro Alligator Ear Forcep
Suction Controller
Rosen Inserter
Grommet

EI0973 Standard

Jobson Horne Probe
Ear Speculas- 4,5,6 MM
Myringotomy Knife
Suction Tube - 18,20 Gauge
Micro Alligator Ear Forcep
Suction Controller
Rosen Inserter
Grommet



Coaglator System



Controller Input

Controller Input

Controller Input

Controller Input

Output Power

Output Power

Controller Dimensions

Controller Dimensions

Controller Dimensions

Controller Dimensions

Voltage

Frequency

RMS Current

Fuse Rating

Fundamental Frequency

Max Output Power

Weight (Max)

Height

Width

Length

220V± 10 %

50Hz ± 1Hz

2 Amp

10A/250V

450 kHz


≤ 300W @ 275Ω

6.5kg

140mm

350mm

350mm



Modus Operands

Working Modes	Functions	Power Levels	
Cutting Mode	Cutting & Vaporization	Power 4-9	6 Adjustable Power Levels
	Cutting & Coagulation	Power 6-8	3 Adjustable Power Levels
	Vaporization & Puncture	Power 3-5	3 Adjustable Power Levels
Coagulation Mode	Ablation & Coagulation	Power 3-5	3 Adjustable Power Levels

Timing Function: Timing for coagulation mode: 0-9 seconds adjustable.

Wand

REF Number

Units/Box, ETO Sterile

Advance Tonsil Wand	EI2302	1 Unit
Advance Adenoid Wand	EI2303	1 Unit
Turbinate Wand (1.3mm)	EI2304	1 Unit
Turbinate Wand (1.6mm)	EI2305	1 Unit
Suction Hemostasis Wand	EI2307	1 Unit
Nasal Hemostasis Wand	EI2306	1 Unit
Larynx Wand	EI2308	1 Unit
Micro Laryngeal Wand	EI2309	1 Unit




Laryngectomy Tube

Material : Medical Grade Silicone

1 Unit/Box
ETO Sterile/Non-Sterile

REF Number Standard	REF Number Fenestrated	Size	Inner Diameter (A) mm	Outer Diameter (B) mm	Length Diameter (L) mm
E3110836S	E3110836F	8/36	9.5	12	36
E3110855S	E3110855F	8/55	9.5	12	55
E3110936S	E3110936F	9/36	10.5	13.5	36
E3110955S	E3110955F	9/55	10.5	13.5	55
E3111036S	E3111036F	10/36	12	15	36
E3111055S	E3111055F	10/55	12	15	55
E3111236S	E3111236F	12/36	13.5	17	36
E3111255S	E3111255F	12/55	13.5	17	55





Laryngectomy Button

Material : Medical Grade Silicone

1 Unit/Box
ETO Sterile/Non-Sterile

REF Number	Size	Inner Diameter (A) mm	Outer Diameter (B) mm
E3131208	12/8	12	8
E3131408	14/8	14	8
E3131608	16/8	16	8
E3131808	18/8	18	8
E3131218	12/18	12	18
E3131418	14/18	14	18
E3131618	16/18	16	18
E3131818	18/18	18	18




Laryngeal Umbrella Keel

Material : Medical grade Silicone

1 Unit/Box
ETO Sterile/Non-Sterile

REF Number Clear	REF Number White	Size (A) mm	Size (B) mm
E31212C	E31212W	12	16
E31214C	E31214W	14	16
E31216C	E31216W	16	16



Heat/Moisture Exchanger (HME)

REF
Number

Descriptions

Units/Box
ETO Sterile/Non Sterile

EA91001	Air-Heat/Moisture Exchanger (HME)	5 Units
EA91002	Moist-Heat/Moisture Exchanger (HME)	5 Units
EA91006	Slim HME	30 Units






Salivary by-pass Tube

Material : Medical grade Silicone

1 Unit/Box
ETO Sterile/Non-Sterile

REF No. Clear	REF NO. White	Size (A) mm	Size (B) mm	Length (C) mm
E31410C	E31410W	10	24	188
E31412C	E31412W	12	24	188
E31414C	E31414W	14	27	188




Standard Tracheal T-Tube

Material : Medical Grade Silicone

1 Unit/Box
ETO Sterile/Non-Sterile

REF No. Clear	REF No. White	Outer Diameter (mm)	A (mm)	B (mm)	C (mm)	D (mm)
E31506CS	E31506WS	6	10	6	48	5
E31508CS	E31508WS	8	16	12	36	6
E31510CS	E31510WS	10	23	17	44	8
E31512CS	E31512WS	12	31	21	55	11
E31514CS	E31514WS	14	32	23	64	11
E31516CS	E31516WS	16	31	23	63	12





Long & Extra Long Tracheal T-Tube

Material : Medical Grade Silicone

1 Unit/Box
ETO Sterile/Non-Sterile

REF No. Clear	REF No. White	Outer Diameter (mm)	A (mm)	B (mm)	C (mm)	D (mm)	Size
E31506CL	E31506WL	6	34	7	48	5	Long
E31508CL	E31508WL	8	58	12	39	7	Long
E31510CL	E31510WL	10	63	17	50	8	Long
E31512CL	E31512WL	12	72	20	59	11	Long
E31514CL	E31514WL	14	72	23	66	11	Long
E31516CL	E31516WL	16	81	24	70	12	Long
E31506CE	E31506WE	6	40	30	48	5	Extra-Long
E31508CE	E31508WE	8	60	49	40	7	Extra-Long
E31510CE	E31510WE	10	68	48	41	8	Extra-Long
E31512CE	E31512WE	12	98	39	64	11	Extra-Long
E31514CE	E31514WE	14	97	48	74	11	Extra-Long



Speakvalve

REF Number	Product Name	Units/Box, Non Sterile
EA91021	Speakvalve	1 Unit



Shower Guard | Laryngeal Brush | Laryngeal Holder

REF Number	Descriptions	Units/Box Sterile/Non-sterile
EA91003	Shower Guard	1 Unit
EA91004	Laryngeal Tube Brush	1 Unit
EA91005	Laryngeal Tube Holder	1 Unit





Absorbable Haemostatic Gelatin Sponge



REF Number	Product	Specification	Units/Box Gamma Sterile
E91018050	Gelatin Sponge	80mm x 50mm x 1mm	20 Units
E91017050	Gelatin Sponge	70mm x 50mm x 1mm	20 Units
E91048050	Gelatin Sponge Bilateral	80mm x 50mm x 4mm	10 Units
E91048025	Gelatin Sponge Unilateral	80mm x 25mm x 4mm	10 Units
E91078020	Gelatin Sponge	60mm x 20mm x 7mm	10 Units
E91078015	Gelatin Sponge Nasal	80mm x 15mm x 7mm	10 Units
E91078025	Gelatin Sponge Gynec	80mm x 25mm x 7mm	10 Units
E91072020	Gelatin Sponge	20mm x 20mm x 7mm	10 Units
E91108050	Gelatin Sponge	80mm x 50mm x 10mm	10 Units
E91107050	Gelatin Sponge	70mm x 50mm x 10mm	10 Units
E91108040	Gelatin Sponge Strifen	80mm x 40mm x 10mm	10 Units
E91108020	Gelatin Sponge Strip	80mm x 20mm x 10mm	10 Units
E9115200	Gelatin Sponge	200mm x 150mm x 15mm	10 Units
E911010	Gelatin Sponge Dental	10mm x 10mm x 10mm	32 Units
E911030	Gelatin Sponge Dial	30mm x 30mm x 10mm	10 Units
E9110120	Gelatin Sponge 120	120mm x 100mm x 10mm	02 Units
E9110125	Gelatin Sponge	125mm x 80mm x 10mm	02 Units
E9110250	Gelatin Sponge 250	250mm x 80mm x 10mm	02 Units
E9105200	Gelatin Sponge Film	200mm x 70mm x 0.5mm	10 Units
E910203	Gelatin Sponge MicroSponge	3mm x 2mm x 2mm	200 Units
E910404	Gelatin Sponge MicroSponge	4mm x 4mm x 4mm	100 Units
E918005	Gelatin Sponge Nasal Tampon	5mm Dia x 80mm	10 Units
E918030	Gelatin Sponge Anal Tampon	30mm Dia x 80mm	06 Units
E911208	Gelatin Sponge Dental Tampon	8mm Dia x 12mm	32 Units
E9101	Gelatin Powder	1 Gram	06 Units

INFRASTRUCTURAL HIGHLIGHTS

50000 sq. feet "state of the art" world class manufacturing infrastructure.

Fully automated production, assembly & packaging lines.

Fully automated Medical Grade Silicone Moulding lines.

Clinically approved biocompatible product portfolio.

Fully automated CNC & injection moulding lines.

1 million units per day production capacity.

In-house Microbiology Culture laboratory.

Class 10000 certified clean room facilities.

Class 1000 Laminar Flow systems.

In-house ETO sterilisation plant.



EON MEDITECH PVT. LTD.

12, Gurukrupa Society, Uttran Road,
Surat - 394105, Gujarat, INDIA.

Ph.: +91 261-2494330

info@eonmeditech.com



EON MEDITECH PVT. LTD.

BL 54, 127, 129 Khanpur Road,
Mirapur, Kamrej, Surat-394320,
Gujarat, India.

Ph.: +91 2621-264330

www.eonmeditech.com



Largest manufacturing infrastructure

1 million units per day production capacity

Presence in 100 countries